




JUST A CLICK

No complicated steps. No unwieldy jargons. Convert all your **DICOM** files into **JPG, BMP** or **AVI** formats in just a click.



Simplicity is the biggest advantage of **Mystiq**. It converts DICOM files into multiple file types including JPG, BMP and AVI. Stores them in predefined directories in an automated manner. Runs as an application service that accepts requests from multiple DICOM sources. And stores them in structured directories as customized by you.



D&Z \ 210108 \ Charles V

HOW Mystiq WORKS

1. Select the desired DICOM images and send it to **Mystiq**
2. As per the user configuration, **Mystiq** will convert the selected DICOM files and store them in specified directory structure.

KEY FEATURES

- One click conversion to JPG, BMP or AVI.
- Easy to install and use.
- Multiple users can concurrently convert DICOM files because the application runs as a service.
- Configurable directory structure to store images based on modality, study date patient ID, accession number, study UID and number of frames (single frame or multi frame).
- Send images from any DICOM source including modality, server, viewer, workstation, etc.
- Supports all DICOM modalities including Ultrasound, CT, MRI, CR, DR, etc.
- Security logs for tracking errors.

SYSTEM REQUIREMENTS

- **Operating System:** Windows 2000 (SP4) or Windows XP (SP2).
- **Processing Power:** Pentium IV 2.0 Ghz or more.
- **Memory:** 512 MB RAM or more.
- **Hard Disk Space:** 50 MB of Free Disk Space for installation.

ASHVA TECHNOLOGIES LTD.

505 Hwy 169 N, Waterford Park, Suite 300, Plymouth, MN 55441.
Ph: 763 746 9223. Fax: 763 746 9224

Ashva Technologies Sdn. Bhd. (742382-M)

No.104 A, Jalan SS
24/2, Taman Megah, 47301 Petaling Jaya, Selangor Darul Ehsan, Malaysia.
H/P: (6017) 2340805

ASHVA TECHNOLOGIES PVT. LTD.

#15 ARK Colony, 'Urmilla House' - 1st Floor, Eldams Road, Alwarpet,
Chennai - 18, India.
Ph: 91-44-4208 8867 / 76. Fax: 91-44-4208 8897



www.ashva.com